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1. ✓ Bclr jALA * 8
- 2.
- 3.
4. ✓
5.
 - a.
 - b.
 - c.


FIG. 19 -

6. Confirm the current message.
- ✓ The machine is ready for operation.

Switching off the machine

Danger of an electric shock if unplugging the power cable before switching off the machine

If you unplug the power cable while the main power switch is still in “ON” position, residual voltage in the power cable may cause you to suffer from an electric shock.

- » **Before** unplugging the power cable, switch off the machine at the main power switch.

To switch off the machine, do the following:

1. Ensure that the working chamber is clean.
2. Switch off the machine at the main power switch.
3. (Optional) Unplug the power cable.
4. (Optional) Switch off the main power switch of your workstation / facility.

Mounting & removing blanks

The machine can process the following blanks:

- Discs with a diameter of 98.5 mm – 98.8 mm
- Blocks, max. size: 45 x 20 x 20 mm (L/D/H)*
- Prefabricated abutments*

*requires extra equipment

- i** You can obtain extra equipment from customer service.

How to mount the different blank types

Blank type	Holder required?	How to mount
Discs	No	Discs >> Blank holder
Blocks	Yes	Blocks >> Block holder >> Blank holder
Prefabricated abutments	Yes	Prefabricated abutments >> Prefabricated abutment holder >> Blank holder

Mounting discs in the working chamber

1. Open the working chamber door.
2. Loosen the 4 screws (marked orange in the figure) which attach the fixing ring to the blank holder.


FIG. 20 – LOOSENING THE SCREWS OF THE FIXING RING

3. Remove the fixing ring by turning it clockwise and pulling it upwards. Remove the blank from the blank holder if any.


FIG. 21 – REMOVING THE FIXING RING FROM THE BLANK HOLDER

4. To insert the disc to be processed, do the following:
 - a. If the blank is a multilayer blank, orientate the blank so that the top color layer is on top.
 - b. If the disc has been processed before, ensure the following:
 - The original top side faces upwards.
 - If the disc has been marked for safer remounting by DentalCNC, ensure that the marking is in the position shown below.
 - c. Put the disc into the blank holder.


FIG. 22 – INSERTING A DISC INTO THE BLANK HOLDER (POSITION FOR SAFE REMOUNTING MARKED BLUE)

5. Place the fixing ring on top of the disc in the blank holder and center it by turning it.
6. Using the delivered torque wrench, screw down the fixing ring.

! Tighten the screws firmly. If blanks move or vibrate during processing, the result may become unusable.


FIG. 23 – ATTACHING THE DISC WITH THE FIXING RING

✓ The disc is mounted in the working chamber

Mounting blocks to the block holder

You can mount up to 3 blocks to the block holder.

- ! The block type, the number of blocks and the block positions must correspond to the corresponding job in DentalCAM.


FIG. 24 – THE K5 BLOCK HOLDER

1. Positioning aids of the block holder
 2. Positioning pins
 3. Holes for block shafts
 4. Fixing screws
1. Use the provided torque wrench to loosen the fixing screw at the desired position.

FIG. 25 – UNSCREWING THE SCREW AT THE DESIRED POSITION

2. If the blank is a multilayer blank, orientate the blank so that the top color layer is on top.
3. Position the block so that the positioning pin at the block holder lies in the groove in the block shaft.


FIG. 26 – POSITIONING PINS (LEFT MARKINGS) MUST LIE IN THE GROOVE IN THE SHAFT (RIGHT MARKING)

4. Insert the shaft of the block into the desired position in the block holder until it is firmly seated.


FIG. 27 – INSERTING A BLOCK INTO THE BLOCK HOLDER

5. Use the provided torque wrench to fixate the block with the corresponding fixing screw. Tighten the screw firmly.

FIG. 28 – SCREWING DOWN THE SCREW TO IMMOBILIZE THE BLOCK

- ✓ You can now mount the block holder.

Mounting the block holder in the working chamber

1. Open the working chamber door.
2. Open the blank holder.
3. Remove the blank from the blank holder if any.
4. Correctly orient the block holder:
 - The fixing screws must face upwards.
 - The positioning aids must be positioned on the left side.
5. Insert the block holder into the blank holder and hold it in place with one hand from below.


FIG. 29 – INSERTING THE BLOCK HOLDER INTO THE BLANK HOLDER

6. Close the blank holder with the other hand.

Using an optional abutment holder

With an optional abutment holder, your machine can process prefabricated abutments with prefabricated connection geometries.

- i** Abutment holders for common abutment systems can be obtained from customer service.

You can find more information about abutment holders and abutment blanks at the following internet address: dentalportal.info/abutments

Managing tools

NOTICE

Damaging of the spindle or the tool positions if you use improper tools

Improper tools can damage the collet chuck of the spindle and / or the tool positions.

- » Only use tools with a sufficiently large chamfer at the tool shank.
- » Install a retaining ring as a stop ring according to DIN 471-A3.
- » Only insert tools with a maximum diameter of 3 mm at the thickest part into the collet chuck.
- » Only insert tools with a maximum cutting edge diameter of 2.6 mm into the tool magazine.

💡 We recommend original tools as they are designed especially for the designated jobs.

You can insert up to 16 tools into the tool magazine.


FIG. 30 – TOOL MAGAZINE IN THE WORKING CHAMBER (MARKED ORANGE)

Normally, the machine loads tools from the tool magazine of the machine. In the following cases, this might not be possible:

- A tool magazine position is worn so that it can no longer properly hold a tool
- You use tools which do not fit into the tool magazine

In these cases, you can force a manual tool change in DentalCNC (see the corresponding documentation).

- i** If tool magazine inserts are worn, you should replace them. ([↗ Replacing the tool magazine inserts – on page 43](#))

Inserting and exchanging tools

You need to insert or exchange tools in the following cases:

- Upon first use of the machine
- After exchanging tool magazine inserts due to wear
- When the tool life of a tool has expired
- When tools are damaged or worn
- If the next jobs require additional / different tools than those in the tool magazine

- i** You can check the remaining tool life of all tools in the **Tool management** view in DentalCNC.

You can equip the tool magazine of your machine in two ways:

- By manually inserting the tools into the tool magazine. This is the fastest way. It is described below.
- Via the spindle – you insert a tool into the collet chuck and the spindle will deposit the tool in the tool magazine. This option takes more time, but might be more convenient for some. It is described in the documentation for the manufacturing software.

You insert tools manually as follows:

1. Remove worn or damaged tools from the tool magazine if any.
2. Insert the new tools:
 - a. Ensure that the positions of the tools in the tool magazine match the tool positions in DentalCNC.
 - b. Insert the tools straight into the tool positions with the cutting edge pointing downwards. Push them in until the ring touches the rubber.

! If the positions of the tools in the tool magazine do not match the tool positions in DentalCNC, the machine will use the wrong tool(s) during job execution and the job result will become unusable.

3. After replacing a tool, reset the tool life value. For this, use the Tool management view in DentalCNC.



FIG. 31 – INSERTING TOOLS INTO THE TOOL POSITIONS

1	U030-R2-40	9	P250-F1-40
2	U050-F2-40	10	Z100-R2-40
3	U060-R2-40	11	Z100-R2D-40
4	U120-F2-40	12	Z120-F2D-40
5	P100-R1-40	13	Z200-R3-40
6	P100-R2-40	14	Z200-R3D-40
7	P200-R1-40	15	Z060-R2D-40
8	P200-R2-40	16	G260-T-35

FIG. 32 – TOP: TOOL POSITIONS 1 – 16 IN THE TOOL MAGAZINE
 BOTTOM: TOOL POSITIONS 1 – 16 IN DENTALCNC

7 Operation: Executing jobs

You control and start manufacturing with DentalCNC. In this chapter, we will give you a brief overview. For the complete instructions, see the documentation for DentalCNC.

NOTICE

Damaging of the machine when using damaged tools or blanks

If tools or blanks are damaged, parts can break off and damage the machine during job execution.

- » Check the blanks and tools **thoroughly** for damage before every job execution.

When you have prepared your jobs and your machine, you can start machining. Machining is a fully automated process and only requires your attention in case of unexpected events.

- ! Do not move the machine during job execution, otherwise the results may become imprecise.

Starting jobs

1. Ensure the following:
 - You created a job on your CAM computer. It has been transferred to DentalCNC.
 - All required tools are in the right positions in the tool magazine and are neither worn nor damaged. They have also been added to the virtual tool magazine in DentalCNC.
 - The required blanks are mounted.
 - The compressed air supply is set correctly.
2. Close the working chamber door.
3. If you manually control the suction device, switch it on and set it to the required level.
4. Start machining via the depicted icon in DentalCNC.

Aborting machining

You can abort machining as follows:

1. Select the depicted icon.
2. Confirm the current message.

- ✓ The following happens:
 - a. Machining stops immediately.
 - b. You are prompted to select whether the tool in the spindle can be inserted safely into the tool magazine.
- 3. If you want the machine to automatically put the tool into the tool magazine, confirm the dialog. If you want to remove the tool manually from the collet chuck, answer the question in the negative.

- i Removing the tool manually is necessary in case that putting it into the tool magazine would damage the tool magazine inserts.

- ✓ Depending on your choice, the spindle either puts the tool into the tool magazine or moves into the tool change position. In the latter case, continue with the next step.

CAUTION

Danger of cuts and burns when touching tools with your bare hands

If you handle tools on the cutting surface, you may be injured. As the tool may be very hot, you may also suffer from skin burns.

- » Only touch tools at their shank.
- » When handling tools, wear protective gloves.

4. Open the working chamber door.
5. Hold the tool in the collet chuck in place.
6. Confirm the current message.

- ✓ The following happens:
 - a. The collet chuck opens.
 - b. The current dialog window closes.
 - c. A dialog window opens.
- 7. Remove the tool from the collet chuck.


FIG. 33 – REMOVING THE TOOL FROM THE COLLET CHUCK

8. Close the working chamber door.
9. Confirm the current message.
- ✓ The following happens:
 - a. The collet chuck closes.
 - b. The working chamber door closes.
 - c. The spindle moves to its default position.

Job interruptions and job abortions

A job will be *interrupted* in the following cases:

- The compressed air supply is not sufficient
- The vacuum in the working chamber is not sufficient

An *interrupted* job will normally be continued automatically after the error is corrected.

A job will be *aborted* in the following cases:

- In case of a machine malfunction
- In case of a tool breakage
- In case of a power failure

If a job was *aborted*, you have to restart it.

i DentalCNC allows you to resume the job from the last machining step. For more information, see the corresponding documentation.

How to proceed in case of a job interruption

If the job was interrupted, DentalCNC displays a corresponding message.

If the compressed air is insufficient

- » Check the following:
- The manometer of the compressed air regulator
 - The installation of the pneumatic hoses
 - Your compressor

If the vacuum is insufficient

- » Check the suction hose and your suction unit.

How to proceed in case of a machine malfunction

A machine malfunction is recognized by the internal control unit in case of a critical event. The working chamber will be illuminated in red. DentalCNC displays the error message and error code that was sent by the control unit.

1. Record the error message and error code that is displayed.
2. Restart the machine and the CAM computer. If the problem persists, continue with the next step.
3. Disconnect the machine from the electrical source and prevent it from being restarted.
4. Contact customer service. Have the error message and error code readily available.
5. If you need to remove a blank from the working chamber, perform an emergency opening of the working chamber door.

How to proceed in case of a tool breakage

If a tool breaks during machining, the machine will not recognize this immediately. Instead, the spindle will continue to move with the broken tool. The tool breakage will be recognized upon the following events:

- The next regular tool change
- The next check for broken tools, if activated in the application settings of DentalCNC.

A tool breakage can be caused by the following:

- The tool was damaged or worn
- The tool was put into the wrong tool magazine position or was manually inserted into the spindle at the wrong time. As a consequence, it was not suitable for the processing step.
- The distribution of the objects in the blank (“nesting”) was not correct.

If a tool breaks, do the following:

1. Open the working chamber door.
2. Remove all parts of the broken tool from the working chamber and the collet chuck.
3. If the spindle picked up the tool from the tool magazine, check if the tool was inserted into the correct position. Insert a spare tool into the correct position in the tool magazine.
4. If you manually inserted the tool into the collet chuck, check if the broken tool corresponds to the

tool type which you were prompted to insert. Have a correct spare tool ready.

5. Close the working chamber door. Restart the job.

 DentalCNC allows you to resume the job from the last machining step. For more information, see the corresponding documentation.

 If tools break regularly, you can find additional information in the troubleshooting section: ([↗ Troubleshooting](#) – on page 48)

How to proceed in case of a power failure

As long as the machine is not powered, you do not have access to the working chamber.

- » After a short power failure, restart the machine and the CAM computer.
- » If you need to access the working chamber in case of a longer power failure, perform an emergency opening of the working chamber door.

Emergency opening of the working chamber door

WARNING

Crushing hazard and cutting injuries if the working chamber door is open

If the working chamber door is open during machining, it will not protect users from bruises and cuts.

- » **Do not** open or close the working chamber door during machining.
- » **Never** operate the machine with the working chamber door open.
- » Perform the emergency opening only if you are authorized to do so and if you have received training on it.
- » Activate the safety interlock **immediately** after completing your work in the working chamber.

CAUTION

Cutting injuries when touching a rotating tool

If a power failure or a machine malfunction occurs during machining, the spindle including the inserted tool keeps rotating. If you touch the rotating tool, you will suffer from cutting injuries.

- » Wait until the spindle has stopped rotating before performing an emergency opening.

The working chamber door is equipped with a safety interlock that only allows you to open the door when the machine is powered. You can use the emergency release of the working chamber door if the working chamber door is permanently locked. There is an opening for the emergency release at the right side of the machine.

FIG. 34 – OPENING FOR THE EMERGENCY RELEASE

You can perform an emergency opening as follows:

1. Switch off the machine at the main power switch. Disconnect the machine from the electrical source.
2. Deactivate the safety interlock of the working chamber door, as shown in the figure (↗ Fig. 35 – on the next page):
 - a. Remove the protective cap from the emergency release opening.
 - b. Insert the key for the emergency release of the working chamber door through the opening.
 - c. Deactivate the safety interlock of the working chamber door by turning the key counter-clockwise by 90°


FIG. 35 – USING THE KEY FOR THE EMERGENCY RELEASE

Top: Remove the cap
Middle: Insert the key
Bottom: Turn the key

3. Open the working chamber door.
4. Carry out your work in the working chamber.
5. Close the working chamber door.
6. Reactivate the safety interlock of the working chamber door by turning the key for the emergency release clockwise by 90°.
7. Check if the working chamber door actually is locked.
8. If you can still open the working chamber door, repeat the 3 previous steps.
9. Close the opening for the emergency release again with the protective cap.

8 Maintenance and do-it-yourself

Some day-to-day basic maintenance and preventive maintenance is essential to keep the machine mechanics and electric components in good condition for proper machining results.

It is your responsibility to make sure that preventive maintenance, as well as basic maintenance, is performed.

You are the only one who can ensure that your machine receives the proper maintenance care. You are a vital link in the maintenance chain.

Maintenance section

For your convenience, DentalCNC lists all basic maintenance tasks in the **Maintenance** section. In the **Maintenance** section, you can see when the individual tasks are due.

When you have marked a maintenance task as complete, its time interval is reset and the list is updated accordingly.

- » When you have performed a maintenance task, mark it as complete in the **Maintenance** section. This will keep the list up-to-date.

Basic maintenance

Basic maintenance includes tasks which are part of everyday operation. You are responsible to ensure that these tasks are carried out according to the maintenance table. You only need minimal manual skills for these tasks and most required tools are provided with the machine.

Preventive maintenance

Preventive maintenance for this machine has to be scheduled every 2 years, or after 2,000 operating hours at the latest.

- » To schedule preventive maintenance, contact customer service.

Where to get service?

Customer service is your main contact for all service related questions. They will provide you with spare parts, maintenance tips and they will perform preventive maintenance for you on request.

- » When your machine is delivered or installed, ask the service technician for contact details for the customer service team. We also recommend scheduling the first preventive maintenance appointment at this point to ensure your machine gets proper maintenance.

Definition of wear parts

The machine and the extra equipment are warranted for a period of 24 months or 2,000 operating hours, whatever comes first. The warranty covers defects in materials or fabrication as long as the regulations for using the machine in all documents are followed.

Of course, the warranty also covers wear parts as long as their failure cannot be attributed to the function-related abrasion. The wear parts that are mentioned in the maintenance table can already wear down within the warranty period due to their normal function. The average useful life of the wear parts can be seen in the maintenance table.

Use these values to determine operating costs, to plan your spare part stock as well as to create individual maintenance and service plans.

Maintenance table

[🔗](#) page 1

Cleaning the working chamber

Cleaning the working chamber includes the following components:

- Measuring key
- Blank holder
- View window

NOTICE

Damaging of the linear guides or the spindle when cleaning with compressed air

If you clean the working chamber with compressed air, material chips can reach the linear guides or the spindle bearings.

- » **Never** clean the working chamber with compressed air.

CAUTION

Breathing difficulties caused by processing dust

Processing dust that gets into your lungs can cause breathing difficulties.

- » Clean the machine only if the air extraction system is properly installed and activated.

Wear a face mask of class FFP2 during the entire cleaning

We recommend cleaning the working chamber after all other required maintenance tasks.

1. Have ready:
 - A dry cloth
 - A vacuum cleaner
 - A wet brush for the measuring key
2. Close the working chamber door.
3. Move the spindle to the cleaning position by selecting the depicted icon in the **Machining** view in DentalCNC. 
4. Open the working chamber door.
5. Suck up the coarse dirt in the working chamber with the vacuum cleaner.
6. Clean all surfaces and cracks in the working chamber thoroughly with the dry cloth. Use a mild cleaning agent if necessary.
7. Clean the measuring key with the wet brush:
 - a. Clean every opening of the protective cage (marked orange) with the wet brush.
 - b. Clean the measuring key from all sides with the wet brush, reaching into the openings of the protective cage.
 - c. Clean the protective cage with a cloth.


FIG. 36 – MEASURING KEY (MARKED BLUE) WITH PROTECTIVE CAGE (MARKED ORANGE)

8. Thoroughly clean the blank holder from all sides with a brush. Especially clean all openings and movable parts of the blank holder.
9. If necessary, clean the view window with a cloth. Use a mild cleaning agent if necessary.
10. Close the working chamber door.
11. Move the spindle to the default position by selecting the depicted icon in the **Machining** view in DentalCNC. 

Cleaning the collet chuck

NOTICE

Damaging of the spindle when cleaning with compressed air

If you clean the collet chuck with compressed air or ultrasound, the spindle bearings can be damaged.

- » Clean the collet chuck **only** with the appropriate service set.


FIG. 37 – SPINDLE SERVICE SET

1. Cleaning brush
2. Knurled nut
3. Tube of collet chuck grease
4. Cleaning cone

To clean the collet chuck do the following:

1. Have the spindle service set ready.
2. Open the working chamber door.
3. Open the collet chuck by selecting the depicted icon in the **Machining** view in DentalCNC. 
4. Put the knurled nut to the spindle with one hand. With the other hand insert the measuring pin into the collet chuck and keep holding it.


FIG. 38 – UNSCREWING THE COLLET CHUCK

5. Loosen the collet chuck with the knurled nut. Unscrew the collet chuck with your hand or with the knurled nut.
6. Remove the measuring pin from the collet chuck. Then put it aside within reach together with the knurled nut.
7. Clean the inner cone of the spindle with the cleaning cone of the service set.


FIG. 39 – CLEANING THE INNER CONE OF THE SPINDLE

8. Clean the collet chuck with the brush of the service set. Make sure that the 3 buffers in the collet chuck do not fall out during cleaning.


FIG. 40 – CLEANING THE COLLET CHUCK

1. Longitudinal slot
2. Buffer
3. Cleaning brush

NOTICE**Damaging of the spindle when using the wrong grease or applying the grease incorrectly**

If you use unsuitable grease or if grease gets into the longitudinal slots of the collet chuck, the machine may get damaged.

- » Ensure that no grease gets into the longitudinal slots of the collet chuck.
 - » Only use a very small, about pinhead-size amount of the grease.
 - » Only use the provided grease of the service set.
9. Put a small amount of the collet grease on the index finger and smear it with the thumb.
 10. Apply the smeared collet grease to the flanks of the collet chuck.


FIG. 41 – GREASING THE COLLET CHUCK; SURFACE IN WHICH GREASE IS TO BE APPLIED IS MARKED IN BLUE; SLOT THAT MAY NOT BE GREASED IS MARKED IN RED

11. Insert the measuring pin into the collet chuck with one hand and keep holding it. Screw the collet

chuck with the knurled nut in your other hand tightly into the spindle.

! Turn the knurled nut as far as you can so that the collet chuck is properly seated in the spindle. Otherwise rotational imperfections can occur during operation which will worsen your processing results

12. Remove the measuring pin from the collet chuck and store it together with the other components of the spindle service set.

Checking the compressed air regulator

NOTICE

Damaging of the machine when compressed air is contaminated

Compressed air that does not fulfill the guidelines for purity according to ISO 8573-1 can damage the machine.

- » Check the water separator of the compressed air regulator daily for contamination.
- » **Never** use the machine if there is water, oil or solid particles in the water separator.

Checking the water separator for condensate

Condensate in the separator usually points to compressed air not being dry enough.

1. Check if water, oil or solid particles piled up in the water separator.
2. If this is the case, switch the machine off immediately and proceed as follows.
3. Disconnect the machine from the compressed air supply.
4. Check the compressed air supply and make sure that the compressed air fulfills the requirements for air purity according to ISO 8573-1. Do not use the machine until the compressed air fulfills this requirement.
5. Drain the water separator by turning the discharge screw counter-clockwise.
- ✓ The condensate is blown out downwards under pressure.
6. Close the discharge screw again by turning it clockwise.


FIG. 42 – OPENING / CLOSING THE DISCHARGE SCREW OF THE COMPRESSED AIR REGULATOR

Exchanging / cleaning the contaminated filter cartridge

You have to clean or exchange the filter cartridge in the water separator in case of strong contamination.

! A strongly contaminated cartridge can lead to a pressure loss.

If the compressed air fulfills the requirements for air purity according to ISO 8573-1, the filter cartridge usually does not have to be changed.

- » If the filter cartridge is contaminated, check the purity of your compressed air.

You exchange or clean the filter cartridge as follows:

1. Disconnect the machine from the compressed air supply.
2. Unscrew the bowl of the water separator.
3. Unscrew the filter screw below the filter cartridge.
4. Pull out the filter cartridge and clean it if necessary.

i A new filter cartridge is available as spare part from customer service.

5. Insert the new or cleaned filter cartridge and reassemble the water separator.

Cleaning the housing

NOTICE

Damaging of the housing when using an unsuitable cleaning agent

If you use an unsuitable cleaning agent and/or cleaning tool for cleaning the machine housing, the surface or the adhesive foil can get damaged.

- » To avoid scratches, only use a microfiber cloth to clean the housing.
- » Take care of the adhesive symbols so that they do not peel off. The adhesive foil is especially sensitive to rubbing and strong cleaning agents.

If the usage of a special cleaning agent is necessary to eliminate certain dirt, we recommend checking the suitability of the cleaning agent at a hidden place of the part first.

1. Clean the surface with a dry microfiber cloth.
2. If some dirt cannot be removed this way, moisten the cloth. Use a pH neutral cleaning agent if necessary.

Exchanging the main fuse

The internal power supply of the machine has a main fuse that is accessible from the outside and can be replaced if necessary.

- » As a replacement fuse, only use a fuse of the following type: T6,3A L250V

i A new main fuse is available as spare part from customer service.

1. Switch off the machine at the main power switch.
2. Remove the power cord from the connection panel.
3. Remove the cover of the fuse.


FIG. 43 – COVER OF THE FUSE (MARKED ORANGE)

4. Remove the defective fuse and replace it with a new fuse.
5. If you do not have a replacement fuse ready, take the replacement fuse from the right side of the fuse cover and put it into the left side of the fuse cover.
6. Remount the fuse cover.

Calibrating the axes

NOTICE

Deterioration of machining results caused by an incorrect calibration

At delivery, your machine is already calibrated. As long as your machining results are accurate, a new calibration is not necessary. A calibration takes much time and will deteriorate the machining results if it is improperly executed.

- » In case of inaccurate machining results, try adjusting the machining conditions first: Check the fixation and quality of the blank and the state of the tool.
- » **Before** calibrating the machine, contact customer service.
- » **Be very careful** when measuring and entering data during calibration. When in doubt, abort the calibration.

By calibrating the machine with test and calibration specimens, it may be possible to improve the machining results.

i The documentation for the manufacturing software contains all information on calibrating the machine. Therefore, you will only find information specific to this machine in this document.

Your machine is delivered with a calibration set. It contains the following parts:

- Calibration blanks which are used to mill calibration and test specimens
- A tool for milling the calibration or test specimens
- A micrometer for measuring the machining precision

You can calibrate the machine as follows:

1. Have the calibration set readily available.
2. Mount the calibration blank into the blank holder.
3. Follow the instructions on calibrating the machine in the documentation for DentalCNC.
4. Keep all parts of the calibration set except used calibration blanks for further use.

Replacing the tool magazine inserts

When tool magazine inserts are worn, they should be replaced. New inserts are delivered without holes for the tools. The holes must be drilled into the inserts with the machine.

- Your machine comes with tool magazine inserts as spare parts and with the drill tool.
- Additional inserts and drill tools are available via customer service.

i The documentation for the manufacturing software contains step-by-step instructions for drilling the holes into the inserts. Below you find a description of how to replace the tool magazine insert in the machine.

You can replace tool magazine inserts as follows:

1. Have the spare tool magazine inserts readily available.
2. Open the working chamber door.
3. Remove all tools from the tool magazine.
4. Unscrew the 4 screws on the upper side of the tool magazine and lift the cover.


FIG. 44 – UNSCREWING THE SCREWS AND LIFTING THE COVER

5. Remove the 2 existing tool magazine inserts and replace them with 2 new ones.


FIG. 45 – NEW TOOL MAGAZINE INSERTS (RIGHT, WITHOUT TOOL POSITIONS) REPLACE OLD INSERTS (LEFT, WITH TOOL POSITIONS)

6. Put the cover back onto the tool magazine and screw it down.
7. Follow the instructions in the documentation for DentalCNC and drill the tool positions using the provided drill bit.

Updating the software and firmware

Updating the manufacturing software and the firmware is an interdependent procedure.

- » Be sure to read all information in this chapter.

NOTICE

Damaging of the control unit when a firmware update is interrupted

If the firmware update is interrupted, the control unit of the machine may become permanently damaged.

- » Only update the firmware if a permanent electrical supply of the machine and the computer is guaranteed.
- » Only update the firmware if the computer's operating system is stable and is free of malware.
- » Only update the firmware if the connection between the computer and the machine is stable. Always use a wired connection during the update.
- » Do not disconnect the machine or the computer from the electrical source nor switch off the machine or computer during a firmware update.
- » Do not close DentalCNC during a firmware update.

Updating the manufacturing software

DentalCAM and DentalCNC are updated regularly.

There are 2 ways to update your existing installation:

- Automatically via the update function in DentalCNC (recommended)
- Manually via an installation package. Ask customer service for the latest installation package.

 The automatic update function is not available for all editions of DentalCAM.

- » To update DentalCAM & DentalCNC, follow the instructions in the documentation for the manufacturing software.

Updating the firmware of the machine

The firmware is the internal control software of your machine. New versions may introduce new functions and improve existing ones. New firmware versions come as part of new versions of DentalCNC.

New firmware versions will be automatically installed after a new DentalCAM & DentalCNC version has been downloaded via the automatic update function. In the following cases, you need to update the firmware manually:

- You do not use the automatic update function.

- No connected machine was found during the automatic download of the new DentalCAM & DentalCNC version.

You can manually update the firmware of your machine as follows:

If DentalCNC contains a new firmware, you will be prompted to update the firmware after DentalCNC connected to the machine or when you start machining. The firmware update should take 5 – 15 minutes to complete.

1. To update the firmware, confirm the message. You can cancel the message to postpone the update, but we recommend updating the firmware immediately.
 - ✓ A window opens.
2. In the current window, select the following icon:
 - 
- ✓ The following happens:
 - a. DentalCNC updates the firmware of the connected machine.
 - b. DentalCNC prompts you to restart the machine.
3. Switch off the machine at the main power switch.
4. Wait 10 seconds.
5. Switch on the machine at the main power switch.
- ✓ The machine references. The machine is ready for operation.

Maintenance table

Several times per day

Task	Recommended interval	Procedure / Parts	Tool illustration
Checking the external compressed air lines for damage	Before switching on the machine	Visual check	
Checking the suction hose for damage	Before switching on the machine	Visual check	

Once per day

Task	Recommended interval	Procedure / Parts	Tool illustration
Checking the compressed air regulator (☞ page 41)	Before work Clean or exchange cartridge in case of visible contamination or every 2 years	If necessary, new filter cartridge	
Daily cleaning of the working chamber (☞ page 38)	After work If soiled	Vacuum cleaner, Brush, Dry cloth	

Once per week

Task	Recommended interval	Procedure / Parts	Tool illustration
Cleaning the collet chuck (☞ page 39)	Once per week In case of rotational imperfections	Spindle service set	

When necessary

Task	Recommended interval	Procedure / Parts	Spare part illustration
Updating the software and firmware (☞ page 44)	When an update is available		
Cleaning the housing (☞ page 41)		Microfiber cloth, Water, Mild cleaning agent (optional)	
Exchanging the main fuse (☞ page 41)		Replacement fuse T6,3A L250V	

Wear parts that you can exchange yourself

Wear part	Recommended interval	Procedure / Parts	Spare part illustration
Collet chuck (☞ page 39)	1,000 operating hours* Every year*	Removing & inserting the collet chucks (old/new) as during cleaning	

Wear part	Recommended interval	Procedure / Parts	Spare part illustration
Tool magazine inserts (↗ page 41)	500 operating hours* Every year*		

Wear parts that customer service exchanges for you

Wear part	Recommended interval	Procedure / Parts	Spare part illustration
Spindle bearings (requires spindle exchange by customer service)	2,000 operating hours*		

*These are recommendation guidelines. Depending on the processing material and how well the machine is cleaned, these values may differ.

9 Disposal

Disposing of the machining residues

When disposing of the machining residues, obey the following regulations.

- » Avoid entry of hazardous machining residues into the soil, water or into sewers.
- » Dispose of the machining residues as described by the manufacturer of the material.
- » Obey the national and local laws of the disposal location in any case.
- » If necessary, have the machining residues disposed of by an approved disposal company.
- » Keep a reference sample of the disposal product for at least 6 months.

Disposing of the machine

The machine must not be disposed of with the residual waste. This is indicated by the icon which depicts a crossed out trashcan. In the European Union (EU), this is in accordance with Directive 2012/19/EU.

We will dispose of the machine at no cost. The owner will bear the costs for disassembly, packaging and transport.

- » Before sending in the machine for disposal, contact customer service.
- » If you dispose the machine yourself, obey the national and local laws of the disposal location.
- » If necessary, have the machine disposed of by an approved disposal company.

Dismantling, transport and packaging

[↗](#) page 9

Copyright

Distribution or duplication of all content is only allowed by written consent of vhf camufacture AG. This includes reproduction by presentation and broadcasting.

This document is published by:

vhf camufacture AG

Lettenstraße 10

72119 Ammerbuch, Germany

dentalportal.info

10 Troubleshooting

In case something is not working as intended, take a look at the following troubleshooting guide.

NOTICE

Machine damage due to improper troubleshooting

In case of improper troubleshooting, your machine may get damaged.

- » If you are unsure of how to perform certain steps during troubleshooting or cannot solve the problems, abort the troubleshooting and contact customer service.

Additional symbols in this chapter

- ❓ Question to narrow down the problem
- 💬 Suggested solution

I cannot open the working chamber door

❓ Is the machine operating?

While the axes are moving, you cannot open the working chamber door.

💬 If applicable:

- » Wait until the machine has finished.

❓ Has a power failure occurred at the installation site of the machine?

💬 If applicable:

- » Depending on the duration of the power failure, restart the machine or perform an emergency opening.

❓ Is electricity available at the installation site of the machine?

💬 If applicable:

1. Connect the machine to the electrical source.
2. Switch on the machine at the main power switch.
3. If the working chamber lighting does not illuminate, check if the power cable is properly seated in the machine and is connected to the electrical source.
4. Try connecting the machine to a different socket.

❓ Is the door blocked?

💬 If applicable:

1. Switch on the machine at the main power switch.
2. If the working chamber lighting does not illuminate, check if the power cable is properly seated in the machine and is connected to the electrical source.
3. If the machine is not powered, perform an emergency opening when necessary.

I have installed all components, started the software but the machine does not reference

❓ Is the working chamber door open?

The machine does not reference with the working chamber door open.

💬 If applicable:

- » Close the working chamber door.

❓ Is an Ethernet cable connected to the machine?

💬 If applicable:

- » Check if the Ethernet cable is properly seated in the connector and is undamaged. If possible, use the provided cable.

❓ Is the working chamber illuminated in red?

In this case a machine malfunction occurred.

💬 If applicable:

1. Restart the machine.
2. If the working chamber continues to be illuminated in red, contact customer service.

My machine doesn't execute any jobs although it is connected to the computer

❓ Is the working chamber door open?

The machine does not reference with the working chamber door open.

💬 If applicable:

- » Close the working chamber door.

The machining results are not satisfactory and / or tools keep breaking

❓ Do the tool positions in the Tools section correspond to the tools in the respective tool magazine?

If not, the machine uses the wrong tools during job execution.

💬 How to check this:

1. In DentalCNC, compare the tool positions in the with the tools in the corresponding tool magazine.
2. Replace the wrong tools in the tool magazine with the correct ones.

❓ Is the blank properly mounted?

💬 How to check this:

- » Remove the blank and remount it. If using blocks: The groove of blocks has to properly sit on the corresponding positioning pin.

? Are bolts, fixing mechanisms, gaps and openings of the blank holder contaminated by processing dust?

☞ If applicable:

- » Clean the mentioned components thoroughly.

? Is the measuring key contaminated?

☞ If applicable:

- » Clean the measuring key with a brush.

? Are the tools worn?

☞ How to check this:

1. Visually inspect all tools.
2. Check the tool life values in DentalCNC.
3. Replace worn tools with new ones.

? Do rings on the tools sit in the groove on the tool shank?

☞ How to check this:

- » Visually inspect all tools and push rings that have moved into the groove again.

? Are the tool magazine inserts worn?

☞ If applicable:

- » Replace the tool magazine inserts with new ones.

? Do the parameters of the job in the software correspond to the parameters of the blank?

☞ How to check this:

- » Make sure that the following parameters of the job and of the blank in the machine correspond to each other. Also make sure that they are suitable for the objects that you want to machine.
- Material type
 - Blank dimensions
 - Indications (types) of the individual objects

? Do you use the latest version of the manufacturing software that is released for the machine?

- » To update DentalCAM & DentalCNC, follow the instructions in the documentation for the manufacturing software.

? Are the object files of sufficient quality?

☞ How to check this:

1. Check the quality of object files (STL files) in your CAD application or an STL viewer. Especially, consider the manufacturer information on wall thickness and border thickness.
2. If necessary, adjust the settings of your scanner and scan application.

? Is the collet chuck of the spindle contaminated or is it loosely seated in the spindle?

☞ If applicable:

1. Clean the collet chuck with the provided spindle service set.
2. When you insert the collet chuck into the spindle, make sure that it is properly seated.

? Did you exchange the collet chuck at the recommended interval?

☞ How to check this:

- » Verify the recommended interval for exchanging the collet chuck in the maintenance table. If necessary, exchange the collet chuck.

The computer interrupts the job and displays that the air pressure is too low

i You can resolve this problem while the job execution is interrupted. DentalCNC continues the job as soon as the problem is resolved.

? Is the compressed air regulator properly set?

- » Set the air pressure at the compressed air regulator to a correct value.

🔗 *Adjusting the air pressure with the compressed air regulator – on page 21*

? Is the error caused by the external compressed air supply system?

☞ How to check this:

1. Close the main external compressed air supply valve.
2. Check if all pneumatic hoses are properly seated in their connections and are undamaged.
3. Check if your compressor is switched on and set correctly.
4. Open all required valves of your compressed air supply system.

? Does the air pressure fluctuate significantly so that jobs are frequently interrupted?

☞ If applicable:

1. Verify that your compressor can *permanently* generate the required air pressure and volume flow.

i Not every compressor is designed for commercial use with dental machines.

2. If necessary, replace your compressor with one that meets the recommendations.

🔗 *Installing the pneumatics – on page 19*

The computer interrupts the job and displays that the vacuum is too low

i You can resolve this problem while the job execution is interrupted. DentalCNC continues the job as soon as the problem is resolved.

? Is the suction unit switched on and operating?

🗨️ How to check this:

1. Check if the suction hose is properly seated in the corresponding opening and is undamaged.
2. If the machine controls the suction unit:
 - a. Check if the switching unit or data cable is properly installed.
 - b. Try operating the suction unit without the switching unit or data cable.
3. Switch on the suction unit.
4. Set the extraction level higher until the machine continues the job.

? Is the filter or container of the suction unit full?

🗨️ If applicable:

- » Insert an empty filter into the suction unit or empty the container.

? Does the suction unit have an automatic tapping function?

🗨️ If applicable:

- » Choose a shorter tapping interval.

I have exchanged the tool magazine inserts and now there are no holes for the tools anymore

Tool magazine inserts are delivered without drilled holes for tools. You drill them with the machine.

- » Use DentalCNC to drill new holes into the new inserts.

Index

A

Abutment holder 30
 Accessories container 11
 Anti-soiling concept 12
 Aspiration 22
 Axes 13

C

CAM computer 12
 Collet chuck 39
 Compressed air
 Insufficient 33
 Compressed air regulator 21, 41

D

Definition of wear parts 37
 Disposal 47
 Machine 47

E

Electric connection 18
 Emergency opening 35

F

FAQ 48
 Firmware update 44
 Frequently Asked Questions 48
 Front side of the machine 10

H

Hose connection 23

I

Identification plate 10, 12
 Implants
 Regulations 8

Installation 15

Air extraction system 22
 CAM computer 24
 Compressed air regulator 21
 Pneumatic hose 20
 Pneumatics 19
 Site 16
 Switching unit 23

J

Job interruption 33

M

Machine
 Malfunction 33
 Start with tool in collet chuck 25
 Machine location 16
 Main fuse 41
 Maintenance 8
 Cleaning the collet chuck 39
 Cleaning the housing 41-42
 Cleaning the working chamber 38
 Compressed air regulator 41
 Exchanging the main fuse 41
 General information 37
 Preventive maintenance 37
 Replacing the tool magazine inserts 43
 Wear parts 37
 Maintenance table 45
 Manufacturing software 8
 Measuring key 38
 Mounting & removing blanks 26
 Multi-machine control 24

O

Operating temperatures 16

P

Pneumatic connection 19-20

Pneumatic hose 19-20
Power connection 18
Power failure 34
Prefabricated abutments
 Regulations 8
Preventive maintenance 37

R

Rear side of the machine 10
Repackaging 9

S

Scope of delivery 15
Sound emission 12
Spindle 8
Starting jobs 32
Storage 9
Suction unit 22
 Requirements 22
Switching unit 23

T

Temperatures 16
Tool breakage 33
Tool magazine inserts 43
Transport 9
Transport lock 18
Troubleshooting 48

U

Unattended operation 8

W

Wear parts 37
Working chamber 11
 Clean 38
Working chamber door 10
 Emergency release 35

Original EC Declaration of Conformity

according to EC directive for machinery 2006/42/EC Annex II A

We,

vhf camufacture AG

Lettenstraße 10
72119 Ammerbuch
Germany

hereby declare expressly that the

Machine: CNC milling machine
Type: K5 Impression ; K5 Impression / Design: K5 X
Serial No.: K5ID200000000 – K5ID399999999 ;
KXID200000000 – KXID399999999

fulfills all the relevant provisions of the following directives:

- **2006/42/EC** Machinery directive
- **2014/30/EU** EMC directive

References of the applied harmonized standards according to article 7 passage 2:

- **EN 614-1:2006 + A1:2009**
- **EN ISO 12100:2010**
- **EN ISO 16090-1:2018**
- **EN ISO 13849-1:2015**
- **EN ISO 13849-2:2012**
- **EN 60204-1:2018**
- **EN 61326-1:2013**
- **EN 61326-2-1:2013**
- **EN IEC 61000-3-2:2019**
- **EN 61000-3-3:2013**

References of other standards:

- **IEC 61010-1:2010 + A1:2016**

The manufacturer undertakes to electronically transmit relevant information on the machinery in response to a reasoned request by the national authorities. Person established in the community who is authorized to compile the technical file:

Dipl.-Ing. (FH) Frank Benzinger
Vorstandsvorsitzender / Chief Executive Officer (CEO)
vhf camufacture AG
Lettenstraße 10
D-72119 Ammerbuch

Ammerbuch, 2020/01/13

(Frank Benzinger, CEO)